




# Fraction Collector EFC-96

For Microdialysis sample



12×8

## Specifications

No.of fractions	: 1 channel 96 specimens
	: 2 channels 48 specimens
	: 4 channels 24 specimens
Sampling time	: 1 sec to 999.59 min (1 sec step)
Display	: 1cons
	: File number currently in use
	: Current status
	: Vial number in use
	: Sampling time and elapsed time
	: Set temperature
	: Sample rack image
Sample holding containers	: 12 rows x 8(arranged X-Y)
Sample holding capacity	: 0.8mL
Cooling mechanism	: Electronic cooling method (Peltier device) 4? to room temperature
Dimensions (WDH)	: 280 x 390 x 320 (mm)
Weight	: Approx. 11.6kg
Power source	: 100 to 240 VAC , 50/60 Hz,10A